## CLEAR-TRAC® COMPLETE
Disposable Cannula System



- Threaded design permits easy insertion and provides excellent retention
- Integrated triple-seal design and suture management system
- Transparent body for easy visualization

## CLEAR-TRAC FLEXIBLE
Disposable Cannula System



- Flexible body material
- Integrated triple-seal design and suture management system
- Capped flow port



# Cannulas, Trocars and Obturators


## Diagnostic Cannula Sets






### Diagnostic Cannula Sets

| Reference # | Description |
|---|---|
| 72200829* | 6.0 mm High Flow Diagnostic Cannula, double-valve, rotatable |
| 72200830 | 6.0 mm High Flow Diagnostic Cannula, single-valve, rotatable |
| 8006653 | O-ring for both cannulae |
| 72201439 | 6.0 mm High Flow Diagnostic Cannula, single-valve, rotatable, for use with Stryker™ scopes |
| 72201440 | 6.0 mm High Flow Diagnostic Cannula, double-valve, rotatable, for use with Stryker scopes |
| 72201627 | 6.0 mm High Flow Diagnostic Cannula, single-valve, rotatable, for use with Storz™ scopes |
| 72201628 | 6.0 mm High Flow Diagnostic Cannula, double-valve, rotatable, for use with Storz scopes |
| 72201693 | 5.0 mm Medium Flow Diagnostic Cannula, single-valve, rotatable |
| 72201712* | 5.0 mm Medium Flow Diagnostic Cannula, double-valve, rotatable |



### 6.0 mm High Flow Diagnostic Cannula Set, double-valve

| Reference # | Description |
|---|---|
| 3868 | 6.0 mm High Flow Diagnostic Cannula Set, double-valve, rotatable |
| **Set includes** | |
| 72200829 | 6.0 mm High Flow Diagnostic Cannula, double-valve, rotatable |
| 4356 | Obturator, conical tip |



### 6.0 mm High Flow Diagnostic Cannula Set, single-valve

| Reference # | Description |
|---|---|
| 7204711 | 6.0 mm High Flow Diagnostic Cannula Set, single-valve, rotatable |
| **Set includes** | |
| 72200830 | 6.0 mm High Flow Diagnostic Cannula, single-valve, rotatable |
| 4356 | Obturator, conical tip |


# Cannulas, Trocars and Obturators

## Operative Cannulas



### Operative Cannulas

| Reference # | Description |
|---|---|
| 3175* | 4.5 mm Operative Cannula, no distal holes, for 4.0 mm arthroscopes with light post on OPPOSITE SIDE of direction of view |
| 3673 | 4.5 mm Operative Cannula, distal holes, for 4.0 mm arthroscopes with light post on OPPOSITE SIDE of direction of view |
| 2704* | 4.5 mm Operative Cannula, distal holes, for 4.0 mm arthroscopes with light post on SAME SIDE of direction of view |
| 3176* | 5.5 mm Operative Cannula, no distal holes, for 4.0 mm arthroscopes with light post on OPPOSITE SIDE of direction of view |
| 7204708 | 5.5 mm Operative Cannula, distal holes, for 4.0 mm arthroscopes with light post on OPPOSITE SIDE of direction of view |
| 2711* | 5.5 mm Operative Cannula, distal holes, for 4.0 mm arthroscopes with light post on SAME SIDE of direction of view |
| ☼ 6900074 | 4.0 mm Operative Cannula for arthroscopes with light post on SAME SIDE of direction of view |



### Operative Cannulas with Fluid Seals

For 4 mm DYONICS◇ Arthroscopes with light post on OPPOSITE SIDE of direction of view, J-Locking interface.

| Reference # | Description |
|---|---|
| 7205366* | 5.5 mm Cannula, outflow valve, distal holes, with fluid seal |
| 7205368 | 5.5 mm Cannula, outflow valve, no distal holes, with fluid seal |
| 7205367 | 5.5 mm Cannula, finger post, distal holes, with fluid seal |
| 7205369 | 5.5 mm Cannula, no distal holes, with fluid seal |



### Operative Cannula Accessories

| Reference # | Description |
|---|---|
| 7204864 | Infusion Adaptor with valve for inflow or outflow. Attaches to cannulas REF 3673, 7204708, 2704, 2711, 3175 and 3176 |
| 4536 | Fluid seals (25 per package) |
| 72200027 | Fluid seals, sterile (25 per package) |

* Item shown

☼ Special, Made to Order products are available to meet the individual needs of surgeons. Please appreciate that longer lead time for delivery should be expected. For more information, please contact your sales representative or Customer Service.

# Cannulas, Trocars and Obturators

## Irrigation Extenders



### Irrigation Extenders

| Reference # | Description |
|---|---|
| 72200027 | Fluid seals |
| 72200831* | High Flow Irrigation Extender, double-valve, rotatable |
| 72200832 | High Flow Irrigation Extender, single-valve, rotatable |
| 72201021 | High Flow Irrigation Extender for cannulas with fluid seals, single-valve, rotatable |



| | |
|---|---|
| 72201022* | High Flow Irrigation Extender for cannulas with fluid seals, double-valve, rotatable |



### Irrigation Extenders

| Reference # | Description |
|---|---|
| 7205370* | Irrigation Extender for cannulas, with fluid seals, with outflow valve |



| | |
|---|---|
| 7205371* | Irrigation Extender for cannulas, with fluid seals, rotatable |

* Item shown

# Cannulas, Trocars and Obturators

## Obturators and Trocars









### Obturators and Trocar

| Reference # | Description |
|---|---|
| 7205376 | Obturator, blunt tip, fits all 4.5 mm, 6.0 mm and 6.0 mm diagnostic cannulas |
| 7204812 | Obturator, blunt tip, fits all 6.0 mm diagnostic cannulas with Storz™ interface |
| 4356 | Obturator, conical tip, fits all 4.5 mm and 5.8 mm diagnostic cannulas |
| 72200402 | Obturator, conical tip, 4.5 mm with 1.5 mm, fits all 4.5 mm and 6.0 mm diagnostic cannulas |
| 7207701 | Obturator, conical tip, fits all 6.0 mm diagnostic cannulas with Stryker™ interface |
| 7204811 | Obturator, conical tip, fits all 6.0 mm diagnostic cannulas with Storz interface |
| 7209965 | Obturator, conical tip, fits all 5.0 mm diagnostic cannulas |
| 1992* | 4.5 mm Obturator, conical tip, fits all 4.5 mm operative cannulas |
| 2209* | 5.5 mm Obturator, conical tip, fits all 5.5 mm operative cannulas |
| 7205783 | 5.5 mm Obturator, conical tip, for cannulas with fluid seals |
| 7205377 | Trocar, sharp tip, fits all 4.5 mm and 6.0 mm diagnostic cannulas |
| 7204810 | Trocar, sharp tip, fits all 6.0 mm diagnostic cannulas with Storz interface |
| 1995* | 4.5 mm Trocar, sharp tip, fits all 4.5 mm operative cannulas |
| 7202210 | 5.5 mm Trocar, sharp tip, fits all 5.5 mm operative cannulas |
| 3801* | Switching Stick, 4.3 mm x 13" |
| ☼ 6900734 | Switching Stick, 2.9 mm (for REF 3781) |
| ☼ 6900295Z | Switching Stick for 3.0 mm |
| 3802 | Switching Sleeve, 4.5 mm x 10" |
| 2787 | Wissinger Rod, 4.5 mm |
| 7205395 | Nitinol Guide Wire, 1.2 mm x 45 cm, (3 per package) |
| ☼ 72202791 | 3.0 mm Luer Cannula |



### Outflow Cannula Set

| Reference # | Description |
|---|---|
| 7204910 | Outflow Cannula Set |
| **Set includes** | |
| 7204863 | 3.0 mm Outflow Cannula with distal holes |
| 7204865 | Obturator, conical tip |
| ☼ 6900365 | 3.0 mm Trocar (SOD) |

* Item shown

☼ Special, Made to Order products are available to meet the individual needs of surgeons. Please appreciate that longer lead time for delivery should be expected. For more information, please contact your sales representative or Customer Service.

# Cannulas, Trocars and Obturators

## Operative Cannula Sets – 4.5 mm



### Light Post – OPPOSITE SIDE

| Reference # | Description |
|---|---|
| 3709 | 4.5 mm Operative Cannula Set for 4.0 mm arthroscopes with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 3673 | 4.5 mm Operative Cannula (2 per set) |
| 1995 | 4.5 mm Trocar, sharp tip |
| 1992 | 4.5 mm Obturator, conical tip |
| 3708 | Irrigation Extender, single-valve, fixed |
| 1986 | Infusion Adaptor |



### Light Post – SAME SIDE

| Reference # | Description |
|---|---|
| 2701 | 4.5 mm Operative Cannula Set for 4.0 mm arthroscopes with light post on SAME SIDE of direction of view |
| **Set includes** | |
| 2704 | 4.5 mm Operative Cannula (2 per set) |
| 1995 | 4.5 mm Trocar, sharp tip |
| 1992 | 4.5 mm Obturator, conical tip |
| 2705 | Irrigation Extender, single-valve, fixed |
| 1986 | Infusion Adaptor |



### Light Post – OPPOSITE SIDE

| Reference # | Description |
|---|---|
| 4188 | 4.5 mm Rotatable Operative Cannula Set for 4.0 mm arthroscopes with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 3673 | 4.5 mm Operative Cannula (2 per set) |
| 1995 | 4.5 mm Trocar, sharp tip |
| 1992 | 4.5 mm Obturator, conical tip |
| 3814 | Irrigation Extender, single-valve, rotatable |
| 1986 | Infusion Adaptor |
| ☼ 6900391Z | 4.5 mm Cannula with 6.0 mm distal holes |

* Item shown

# Cannulas, Trocars and Obturators

## Operative Cannula Sets – 5.5 mm



### Light Post – SAME SIDE

| Reference # | Description |
|---|---|
| 3112 | 5.5 mm Operative Cannula Set for 4.0 mm arthroscopes with light post on SAME SIDE of direction of view |
| **Set includes** | |
| 2711 | 5.5 mm Operative Cannula (2 per set) |
| 7202210 | 5.5 mm Trocar, sharp tip |
| 2209 | 5.5 mm Obturator, conical tip |
| 1986 | Infusion Adaptor |



### Light Post – OPPOSITE SIDE

| Reference # | Description |
|---|---|
| 7204709 | 5.5 mm Operative Cannula Set for 4.0 mm arthroscopes with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 7204708 | 5.5 mm Operative Cannula (2 per set) |
| 7202210 | 5.5 mm Trocar, sharp tip |
| 2209 | 5.5 mm Obturator, conical tip |
| 1986 | Infusion Adaptor |



### Infusion Adaptor

| Reference # | Description |
|---|---|
| 1986 | Infusion Adaptor |
| 72200623 | Infusion Adaptor with FMS inflow port (SOD) |

# Hip Arthroscopy Set
## Cannulas and Obturators



ARTHROGARDE® Hip Access Cannulas are designed to be atraumatic, featuring smooth ends for optimal safety and high-strength design for enhanced durability.

The uniquely designed rounded ends provide a smooth surface to help protect soft tissue when performing hip arthroscopy. Plus, they're extra strong to prevent bending or wear and tear over time. Combine with our CROSSTRAC® Hip Guide System for accurate trajectories and first-time access to the hip joint.

Please note, all ARTHROGARDE Cannulas utilize J-SNAP Irrigation Extenders 72202448 or 72202449 at top of next page.

### ARTHROGARDE Hip Access Cannulas

| Reference # | Description |
|---|---|
| 72201741 | 4.5 mm ARTHROGARDE Hip Access Cannula, green (CROSSTRAC Hub) |
| 72202003 | 4.5 mm ARTHROGARDE Hip Access Cannula, green (Standard Hub) |
| 72201742 | 5.0 mm ARTHROGARDE Hip Access Cannula, blue (CROSSTRAC Hub) |
| 72202004 | 5.0 mm ARTHROGARDE Hip Access Cannula, blue (Standard Hub) |
| 72201743 | 5.5 mm ARTHROGARDE Hip Access Cannula, red |
| 7209345 | 3.0 mm Cannula, inflow |
| 72203133 | ARTHROGARDE Slotted Cannula (4.4" working length) |
| 72201590 | ARTHROGARDE Slotted Cannula (6.2" working length) |
| ☼ 72203327 | 5.0 mm Extended Length ARTHROGARDE Hip Access Cannula, blue (Standard Hub) |
| ☼ 72203326 | 5.5 mm Extended Length ARTHROGARDE Hip Access Cannula, red (Standard Hub) |



### Arthroscopic Hip Set Instrumentation – Obturators

| Reference # | Description |
|---|---|
| 72200884 | 4.5 mm Obturator, conical, green, used with REF 72201741 and 72202003 |
| 72200885 | 4.5 mm Obturator, cannulated, green, used with REF 72201741 and 72202003 |
| 72200096 | 4.5 mm Obturator, cannulated, green, large Delrin handle, used with REF 72201741 and 72202003 |
| 72200882 | 5.0 mm Obturator, conical, blue, used with REF 72201742 and 72202004 |
| 72200883 | 5.0 mm Obturator, cannulated, blue, used with REF 72201742 and 72202004 |
| 72200097 | 5.0 mm Obturator, cannulated, blue, large Delrin handle, used with REF 72201742 and 72202004 |
| ☼ 72203329 | 5 mm Extended Length Obturator, cannulated, blue, large Delrin handle, used with REF 72203327 |
| 72200880 | 5.5 mm Obturator, conical, red, used with REF 72201743 |
| 72200881 | 5.5 mm Obturator, cannulated, red, used with REF 72201743 |
| 72200098 | 5.5 mm Obturator, cannulated, red, large Delrin handle, used with REF 72201743 |
| ☼ 72203328 | 5.5 mm Extended Length Obturator, cannulated, red, large Delrin handle, used with REF 72203326 |
| 7207751 | 5.5 mm Obturator, long, cannulated, used with disposable cannula REF 7205360 |
| 7210392 | 5.5 mm Obturator, cannulated, used with slotted cannula REF 72201590 |
| 7209349 | Obturator for 3.0 mm cannula |
| 72201355 | Reuseable Obturator for 5.7 mm Flexible Cannula (use with PN 7205360, 5.7 mm cannula) |
| 7207934 | Hip Arthroscopy System Sterilization Tray (21" x 10" x 3.5") |

# Hip Arthroscopy Set
## Irrigation Extender



Arthroscopic Hip Set Instrumentation – Irrigation Extender

| Reference # | Description |
|---|---|
| 72202448 | J-SNAP quick release Irrigation Extenders, single-valve |
| 72202449 | J-SNAP quick release Irrigation Extenders, double-valve |
| 72200027 | Fluid Seals, sterile (25 per package) |

# Extremity Cannula Sets

## Extremity Cannula Sets



### 2.2 mm Operative Cannula Set for 1.9 mm SHORT SCOPE

| Reference # | Description |
|---|---|
| 4294 | 2.2 mm Operative Cannula Set for 1.9 mm SHORT SCOPE with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 3865 | 2.2 mm Short Cannula, with flow port |
| 3866 | 2.2 mm Short Trocar, sharp tip |
| 3867 | 2.2 mm Short Obturator, conical tip |
| 3312 | TMJ/Wrist Probe |



### 2.9 mm Operative Cannula Set for 2.7 mm LONG SCOPE

| Reference # | Description |
|---|---|
| 3713 | 2.9 mm Operative Cannula Set for 2.7 mm LONG SCOPE with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 3672 | 2.9 mm Short Cannula with finger post (2 per set) |
| 3785 | 2.9 mm Short Trocar, sharp tip |
| 3807 | 2.9 mm Short Obturator, conical tip |
| 3671 | Irrigation Extender, single-valve, fixed |
| 3696 | Infusion Adaptor |
| 72203362 | 2.9 mm Blunt Obturator |




### 2.9 mm Operative Cannula Set for 2.7 mm SHORT SCOPE

| Reference # | Description |
|---|---|
| 3800 | 2.9 mm Operative Cannula Set for 2.7 mm SHORT SCOPE with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 3781 | 2.9 mm Short Cannula with flow port (2 per set) |
| 3807 | 2.9 mm Short Obturator, conical tip |
| 3786 | 2.9 mm Short Obturator, blunt tip |
| 3312 | TMJ/Wrist Probe |

# Extremity Cannula Sets

## Extremity Cannula Sets



### 2.9 mm Operative Cannula Set for 2.7 mm LONG SCOPE

| Reference # | Description |
|---|---|
| 3714 | 2.9 mm Diagnostic Cannula Set for 2.7 mm LONG SCOPE with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 4419 | 2.9 mm Cannula with flow port |
| 3691 | 2.9 mm Trocar, sharp tip |
| 3692 | 2.9 mm Obturator, conical tip |
| 72200411 | 2.9 mm Cannula, single-valve |



### 3.6 mm Operative Cannula Set for 2.7 mm SHORT SCOPE

| Reference # | Description |
|---|---|
| 4317 | 3.6 mm Operative Cannula Set for 2.7 mm SHORT SCOPE with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 4308 | 3.6 mm Operative Cannula, single-valve, fixed |
| 4309 | 3.6 mm Sharp Trocar |
| 4310 | 3.6 mm Obturator, conical tip (accomodates switching stick 6900734) |


Special, Made to Order products are available to meet the individual needs of surgeons. Please appreciate that longer lead time for delivery should be expected. For more information, please contact your sales representative or Customer Service.

# Reusable Cannulas




## Reusable Cannulas

| Reference # | Description |
|---|---|
| 4482 | 4.0 mm Arthroscope Inflow Cannula Set for DYONICS® Arthroscopes with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 4355 | 6.0 mm Rotatable Cannula with finger post |
| 4356 | Obturator, conical tip |



| Reference # | Description |
|---|---|
| 4549 | 4.0 mm Arthroscope Inflow Cannula Set with flush valve for DYONICS Arthroscopes with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 4537 | 6.0 mm Rotatable Cannula with flush valve |
| 4356 | Obturator, conical tip |



| Reference # | Description |
|---|---|
| 4481 | 2.7 mm Arthroscope Inflow Cannula Set standard length for DYONICS Arthroscopes with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 4438 | 4.6 mm Rotatable Cannula with finger post |
| 4439 | Obturator, conical tip |

| Reference # | Description |
|---|---|
| ☼ 6900392 | 2.9 mm Rotatable Cannula, low flow, double-valve (use with 2.7 mm standard arthroscope) |
| ☼ 6900284 | 2.9 mm Obturator for 2.9 mm Cannula (6900392) |

## Diagnostic Cannulas, Extended Length

For use with hip arthroscopy scopes REF 7211054 and 7211055


| Reference # | Description |
|---|---|
| 7211094 | 5.8 mm diagnostic cannula, with inflow and valve, rotatable |
| 7207898 | 4.5 mm Cannulated obturator for cannulas 7211075 and 7207913 |
| 8011677 | FMS port O-ring |
| ☼ 7207913 | 5.8 mm Diagnostic cannula, double-valve, rotatable |

# Reusable Cannulas



## Reusable Cannulas

| Reference # | Description |
|---|---|
| 4525 | 2.7 mm Arthroscope Inflow Cannula Set, short. For DYONICS® Arthroscopes with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 4526 | 3.8 mm Rotatable Cannula with finger post |
| 4527 | Obturator, conical tip |



| Reference # | Description |
|---|---|
| ☼ 7207789 | 3.8 mm Rotatable Cannula, single-valve, short |
| ☼ 7205984 | 3.8 mm Rotatable Cannula with inflow and valve |
| ☼ 7210715 | 3.8 mm Rotatable Cannula, double-valve, short |

| Reference # | Description |
|---|---|
| 7205176 | 4.0 mm Storz™/Zimmer™/Linvatec™ Arthroscope Inflow Cannula Set for 4 mm arthroscopes with light post on OPPOSITE SIDE of direction of view |
| **Set includes** | |
| 7205874 | 6.0 mm Rotatable Cannula with inflow spring lock interface and flush valve |
| 7204810 | Trocar |
| 7204811 | Obturator, conical tip |
| 7204812 | Obturator, blunt tip |



| Reference # | Description |
|---|---|
| 7204479 | Non-scope Inflow Cannula |
| **Set includes** | |
| 4359 | 3.0 mm Inflow Cannula |
| 7204478 | Obturator, conical tip |

## Inflow Cap for Cannulas

| Reference # | Description |
|---|---|
| 7209351 | Inflow cap for cannulas, metal |
| 72202442 | Portal cap, Luer lock interface |
| ☼ 72202791 | 3.0 mm Luer Cannula |

☼ Special, Made to Order products are available to meet the individual needs of surgeons. Please appreciate that longer lead time for delivery should be expected. For more information, please contact your sales representative or Customer Service.

# Cannulas for Use with Arthroscopes

## Use with 2.7 mm Arthroscopes, standard length

| Reference # | Description |
|---|---|
| 72200121 | 4.6 mm Rotatable Cannula, single-valve, standard length |
| 7205999 | 4.6 mm Rotatable Cannula, double-valve, standard length |
| 4439 | Obturator, Conical Tip for 4.6 mm Rotatable Cannulas, standard length |
| 72203536 | 3.8 mm Rotatable Cannula, single-valve, standard length |
| 72201792 | 3.8 mm Rotatable Cannula, double-valve, standard length |
| 7205993 | Obturator, Conical Tip for 3.8 mm Rotatable Cannulas, standard length |
| 72203624 | 2.9 mm Rotatable Cannula, double-valve, standard length |
| 72203625 | Obturator, Conical Tip for 2.9 mm Rotatable Cannulas, standard length |

## Use with 2.7 mm Arthroscopes, short length

| Reference # | Description |
|---|---|
| 7207789 | 3.8 mm Rotatable Cannula, single-valve, short length |
| 7210715 | 3.8 mm Rotatable Cannula, double-valve, short length |
| 4527 | Obturator, Conical Tip for 3.8 mm Rotatable Cannulas, short length |
| 6900513 | Obturator, Blunt Tip for 3.8 mm Rotatable Cannulas, short length |
| 72203138 | Trocar, Sharp Tip for 3.8 mm Rotatable Cannulas, short length |
| 72201751 | 2.9 mm Rotatable Cannula, double-valve, short length |
| 72201752 | Obturator, Conical Tip for 2.9 mm Rotatable Cannulas, short length |

## Use with 1.9 mm Arthroscopes, short length

| Reference # | Description |
|---|---|
| 72201849 | 2.2 mm Rotatable Cannula, double-valve, short length |
| 72201197 | Obturator, Conical Tip for 2.2 mm Rotatable Cannulas, short length |

# Access Disposable Cannulas

## CLEAR-TRAC◇ COMPLETE Disposable Cannula System



### CLEAR-TRAC COMPLETE Disposable Cannula System

The CLEAR-TRAC COMPLETE Cannula System is enhanced with a robust triple-seal design and an integrated suture management system. These new features, combined with the many existing facets from the previous generation, provide a truly comprehensive cannula solution.

- Integrated triple-seal design
- Suture management system
- Removable cannula cap
- Tall thread configuration
- Transparent polycarbonate body
- Latex-free design

### CLEAR-TRAC COMPLETE Cannula Family (sterile, single-use, 10 per box)

| Reference # | Description |
|---|---|
| 72200902 | 8.5 mm x 90 mm Threaded Cannula, with disposable obturator, green, sterile |
| 72200903 | 8.5 mm x 72 mm Threaded Cannula, with disposable obturator, green, sterile |
| 72200904 | 8.5 mm x 45 mm Threaded Cannula, with disposable obturator, green, sterile |
| 72200905 | 7.0 mm x 72 mm Threaded Cannula, with disposable obturator, grey, sterile |
| 72200906 | 7.0 mm x 72 mm Smooth Cannula, with disposable obturator, grey, sterile |
| 72200907 | 5.5 mm x 72 mm Threaded Cannula, with disposable obturator, blue, sterile |
| 72200908 | 5.5 mm x 45 mm Threaded Cannula, with disposable obturator, blue, sterile |
| 72200909 | 4.5 mm x 72 mm Threaded Cannula, with disposable obturator, purple, sterile |



### CLEAR-TRAC COMPLETE Reusable Obturators

| Reference # | Description |
|---|---|
| 72200910 | 8.5 mm x 90 mm Reusable Obturator (4.3 mm cannulation), green |
| 72200913 | 8.5 mm x 90 mm Reusable Solid Obturator, green |
| 72200911 | 8.5 mm x 72 mm Reusable Obturator (4.3 mm cannulation), green |
| 72200912 | 8.5 mm x 72 mm Reusable Obturator (1.5 mm cannulation), green |
| 72200914 | 8.5 mm x 72 mm Reusable Solid Obturator, green |
| 72200915 | 7.0 mm x 72 mm Reusable Obturator (4.3 mm cannulation), grey |
| 72200916 | 7.0 mm x 72 mm Reusable Obturator (1.5 mm cannulation), grey |
| 72201796 | 5.5 mm Reusable Obturator (4.3 mm cannulation), blue |



**Accessories**

| | |
|---|---|
| 3801 | Switching Stick, 4.3 mm, not cannulated (13" long) |

* Item shown

☼ Special, Made to Order products are available to meet the individual needs of surgeons. Please appreciate that longer lead time for delivery should be expected. For more information, please contact your sales representative or Customer Service.

# Access Disposable Cannulas

## CLEAR-TRAC◇ FLEXIBLE Disposable Cannula System



### CLEAR-TRAC FLEXIBLE Disposable Cannula System

The CLEAR-TRAC FLEXIBLE Cannula System is made from a flexible body material, it provides excellent tactile feedback to help identify structures and may reduce the risk of significant tissue damage. Enhanced with a robust triple-seal design. CLEAR-TRAC FLEXIBLE cannulas also feature an integrated suture management system.

- Integrated triple-seal design
- Suture management system
- Flexible body material
- Easy to thread, secure configuration
- Transparent body
- Capped flow port with optional disposable stopcock
- Latex-free design

### CLEAR-TRAC FLEXIBLE Cannula Family (sterile, single-use, 10 per box)

| Reference # | Description |
|---|---|
| 72200424 | 8.0 mm x 90 mm Threaded Cannula, with disposable obturator, lime green, sterile |
| 72200425 | 8.0 mm x 72 mm Threaded Cannula, with disposable obturator, lime green, sterile |
| 72200426 | 8.0 mm x 72 mm Smooth Cannula, with disposable obturator, lime green, sterile |
| 72200427 | 6.5 mm x 72 mm Threaded Cannula, with disposable obturator, orange, sterile |
| 72200428 | 6.5 mm x 72 mm Smooth Cannula, with disposable obturator, orange, sterile |
| 72200429 | 5.5 mm x 72 mm Threaded Cannula, with disposable obturator and trocar, blue, sterile |
| 72200430 | 5.5 mm x 55 mm Threaded Cannula, with disposable obturator and trocar, blue, sterile |
| 72200431 | 5.5 mm x 55 mm Smooth Cannula, with disposable obturator and trocar, blue, sterile |
| 72200432 | 5.5 mm x 45 mm Threaded Cannula, with disposable obturator and trocar, blue, sterile |
| 72200433 | 4.5 mm x 72 mm Threaded Cannula, with disposable obturator, purple, sterile |



### CLEAR-TRAC FLEXIBLE Reusable Obturators

| Reference # | Description |
|---|---|
| 72200899 | 8.0 mm x 90 mm Reusable Obturator (4.3 mm cannulation), lime green |
| 72200900 | 8.0 mm x 72 mm Reusable Obturator (4.3 mm cannulation), lime green |
| 72201796 | 5.5 mm Reusable Obturator (4.3 mm cannulation), blue |
| 72201846 | 6.5 mm x 72 mm Reusable Obturator (4.3 mm cannulation) |



**Accessories**

| | |
|---|---|
| 72200658 | CLEAR-TRAC Disposable Stopcock |
| 3801 | Switching Stick, 4.3 mm, not cannulated (13" long) |

* Item shown

† Disposable knives are supplied as a box of 5, single-use, sterile packs.

# Disposable Cannulas

## CLEAR-TRAC◇ Cannulas



### CLEAR-TRAC Scope/Shaver Cannula

The CLEAR-TRAC Disposable Scope/Shaver Cannula has been specifically designed with users of standard metal operative cannulas in mind. The unique hub design allows irrigation extenders and shaver blades to be connected directly to the cannula, while the integrated double-seal system keeps you dry.

| Reference # | Description |
|---|---|
| 72201634 | CLEAR-TRAC Disposable Scope Shaver Cannula, smooth |
| 72201658 | CLEAR-TRAC Disposable Scope Shaver Cannula, with rib |



### CLEAR-TRAC Screw Cannula

Provides graft protection and aids screw delivery during graft fixation.

| Reference # | Description |
|---|---|
| 7210118* | 7.0 mm CLEAR-TRAC Screw Cannula, latex seal |
| 7210119 | 8.0 mm CLEAR-TRAC Screw Cannula, latex seal |
| 7210120 | 9.0 mm CLEAR-TRAC Screw Cannula, latex seal |
| 7207227 | 7.0 mm CLEAR-TRAC Screw Cannula, latex-free seal |
| 7207228 | 8.0 mm CLEAR-TRAC Screw Cannula, latex-free seal |
| 7207229 | 9.0 mm CLEAR-TRAC Screw Cannula, latex-free seal |

*Single-use, sterile (10 per box)*



### CLEAR-TRAC Threaded Cannula with Obturator

| Reference # | Description |
|---|---|
| 014718* | CLEAR-TRAC Threaded Cannula, 8.0 mm I.D., 76 mm long, green, latex-free seal |
| 014719* | CLEAR-TRAC Threaded Cannula, 5.0 mm I.D., 76 mm long, blue, latex-free seal |
| 7210121 | CLEAR-TRAC Threaded Cannula, 8.0 mm I.D., 76 mm long, green, latex seal |
| 7210122 | CLEAR-TRAC Threaded Cannula, 5.0 mm I.D., 76 mm long, blue, latex seal |

*Includes obturator. Single-use sterile packs (10 per box)*



### CLEAR-TRAC Threaded Cannula with Cannulated Obturator

| Reference # | Description |
|---|---|
| 7207197 | CLEAR-TRAC Threaded Cannula, 10 mm I.D., 76 mm long with 4.0 mm cannulated obturator, white, contains latex |

*Includes cannulated obturator. Single-use sterile packs (10 per box)*

* Item shown

# Disposable
## Cannulas



### Disposable Cannula Sets with Distal Holes

| Reference # | Description |
|---|---|
| 4615 | 4.5 mm Cannula and Obturator, conical tip, distal holes, blue |
| 4616 | 5.5 mm Cannula and Obturator, conical tip, distal holes, gray |
| 4617 | 7.0 mm Cannula and Obturator, conical tip, distal holes, forest green |

*70 mm working length, threaded with seals, latex free (6 per box)*



### Disposable Cannula Sets without Distal Holes

| Reference # | Description |
|---|---|
| 7204895 | 5.5 mm Cannula and Obturator, conical tip, no holes, orange |
| 7204896 | 7.0 mm Cannula and Obturator, conical tip, no holes, lime green |

*70 mm working length, threaded with seals, latex free (6 per box)*



### Universal Cannula – Includes Obturator and Trocar

| Reference # | Description |
|---|---|
| 014623 | 5.0 mm I.D., 55 mm long, blue, latex-free seal |
| 012405 | 5.0 mm I.D., 76 mm long, blue, latex-free seal |
| 012421 | 7.0 mm I.D., 76 mm long (includes obturator only), green, latex-free seal |
| 7210126 | 5.0 mm I.D., 76 mm long, blue, latex seal |

*Single-use sterile packs (10 per box)*

# Cannulae

## 5.5 mm Cannulae, color-coded blue, solid body, no fluid seals

**(same style as fluid seal cannulae but no seal needed)**

| Reference # | Description |
|---|---|
| ☼ 7209953 | 5.5 mm Cannula with distal holes and finger post same side as direction of view |
| ☼ 7209954 | 5.5 mm Cannula with finger post same side as direction of view (no holes) |
| ☼ 7209957 | 5.5 mm Cannula with distal holes and finger post opposite direction of view |
| ☼ 7209958 | 5.5 mm Cannula with finger post opposite direction of view (no holes) |
| ☼ 7210384 | 5.5 mm Blunt Obturator for 5.5 mm color-coded blue cannulae |
| ☼ 72200013 | 5.5 mm Conical Obturator with threaded interface for 5.5 mm color-coded blue cannula |
| ☼ 72200014 | 5.5 mm Trocar with threaded interface for 5.5 mm color-coded blue cannula |
| ☼ 72200012 | 5.5 mm Blunt Obturator with threaded interface for 5.5 mm color-coded blue cannula |

## 5.0 mm Cannulae, color-coded red, solid body, no fluid seals

**(same style as fluid seal cannulae but no seal needed)**

| Reference # | Description |
|---|---|
| ☼ 7209908 | 5.0 mm Cannula with finger post same side as direction of view (no holes) |
| ☼ 7209956 | 5.0 mm Cannula with distal holes and finger post opposite direction of view |
| ☼ 7209857 | 5.0 mm Cannula with finger post opposite direction of view (no holes) |
| ☼ 7209923 | 5.0 mm Trocar for 5.0 mm color-coded red cannulae |
| ☼ 7210385 | 5.0 mm Blunt Obturator for 5.0 mm color-coded red cannulae |
| ☼ 72200010 | 5.0 mm Conical Obturator with threaded interface for 5.0 mm color-coded red cannula |
| ☼ 72200011 | 5.0 mm Trocar with threaded interface for 5.0 mm color-coded red cannula |
| ☼ 72200009 | 5.0 mm Blunt Obturator with threaded interface for 5.0 mm color-coded red cannula |

## 2.7 mm Arthroscope Cannulae (standard length)

| Reference # | Description |
|---|---|
| ☼ 7205989 | 3.8 mm Cannula, FMS, rotatable, with flush-valve |
| ☼ 6900154 | 4.6 mm Cannula, FMS, rotatable, with flush-valve. Used with Obturator (4439) |

## 2.7 mm Arthroscope Cannulae (short length)

| Reference # | Description |
|---|---|
| ☼ 7205984 | 3.8 mm Cannula, FMS, rotatable, with flush-valve. Used with Obturator (4527) |

☼ Special, Made to Order products are available to meet the individual needs of surgeons. Please appreciate that longer lead time for delivery should be expected. For more information, please contact your sales representative or Customer Service.

# Cannulae

## Cannulae

| Reference # | Description |
|---|---|
| 6900022Z | 2.9 mm Cannula, short, with single-valve |
| 4491 | 4.5 mm Diagnostic Cannula, single-valve, rotatable |
| 7205372 | 4.5 mm Diagnostic Cannula, double-valve, fixed body |
| 3134300 | 4.5 mm, no holes, finger post opposite direction of view |
| 6900215Z | 4.5 mm Cannula, operative, with finger post same side as direction of view |
| 6900237Z | 5.0 mm Cannula with distal holes and finger post opposite direction of view |
| 6900084Z | 5.5 mm Cannula with distal holes and finger post same side as direction of view |
| 6900167Z | 5.5 mm Cannula with distal holes and finger post opposite direction of view |
| 6900092Z | 5.5 mm Cannula with finger post opposite direction of view (no holes) |
| 6900106Z | 5.5 mm Cannula with finger post same side as direction of view (no holes) |
| 72201374 | 6.0 mm High Flow Cannula, diagnostic, rotatable, single-valve, with Storz™ spring loaded interface |
| 72201375 | 6.0 mm High Flow Cannula, diagnostic, rotatable, double-valve, with Storz spring loaded interface |

## Trocar

| Reference # | Description |
|---|---|
| 6900279Z | 5.5 mm Trocar for use with 5.5 mm Cannulae |

## Irrigation Extenders

| Reference # | Description |
|---|---|
| 7209959 | Irrigation Extender for color-coded cannula sets, single valve, rotatable |
| 4493 | Irrigation Extender, double-valve, rotatable |
| 72201451 | Irrigation Extender, single-valve, rotatable, removable |

## Miscellaneous

| Reference # | Description |
|---|---|
| 6900274 | 3.0 mm Inflow Cannula with six distal holes. Used with Obturator (7204478) |
| 7210284 | 3.0 mm Inflow Cannula, 90°, no obturator |

# Replacement O-Rings
## Quick Reference Chart

| Cannulas | O-Ring: Scope Side | O-Ring: Cannula Side | O-Ring: FMS Port |
|---|---|---|---|
| 3671 | 8006517 | | |
| 3672 | 8006517 | | |
| 3781 | 8006517 | | |
| 3865 | 8010199 (2 each) | | |
| 3869 | 8006653 | | |
| 4308 | 8003734 | | |
| 4355 | 8006653 | | 8011677 |
| 4359 | | | 8011677 |
| 4419 | 8006517 | | |
| 4438 | 8010199 | | 8011677 |
| 4491 | 8006653 | | |
| 4492 | 8006653 | | |
| 4526 | 8010199 | | 8011677 |
| 4537 | 8006653 | | 8011677 |
| 7205372 | 8006653 | | |
| 7205874 | | | 8011677 |
| 72200829 | 8006653 | | |
| 72200830 | 8006653 | | |

| Cannulas | O-Ring: Scope Side | O-Ring: Cannula Side | O-Ring FMS Port |
|---|---|---|---|
| 2705 | 8006653 | 8006652 (2 each) | |
| 3708 | 8006653 | 8006652 (2 each) | |
| 3814 | 8006653 | 8006652 (2 each) | |
| 4493 | 8006653 | 8006652 (2 each) | |
| 7205370 | 8006653 | 8006652 (2 each) | |
| 7205371 | 8006653 | 8006652 (2 each) | 8011677 |
| 4355 | 8006653 | 8006652 (2 each) | 8011677 |
| 72200083 | 8006653 | 8006652 (2 each) | |
| 72200084 | 8006653 | 8006652 (2 each) | |
| 72200004 | 8006653 | 8006652 (2 each) | 8011677 |
| 4491 | 8006653 | | |
| 4492 | 8006653 | | |
| 4526 | 8010199 | | 8011677 |
| 4537 | 8006653 | | 8011677 |
| 7205372 | 8006653 | | |
| 7205874 | | | 8011677 |
| 72200829 | 8006653 | | |
| 72200830 | 8006653 | | |

# ARTHROPAK
## Customized Procedure Kit




ARTHROPAK Customized Procedure Kit

Smith & Nephew DYONICS® ARTHROPAK Customized Procedure Kits provide increased operational efficiencies and reduced supply expense. You customize your kit with your selection of disposable Smith & Nephew products for any arthroscopic procedure at a reduced cost.

- Free customized service
- Reduced product cost
- Save valuable time and demand on staff
- Streamlined ordering process
- Increased standardization
- Reduced inventory

DYONICS ARTHROPAK Customized Procedure Kits reduce the number of SKUs that need to be ordered and aggregate purchasing volume of arthroscopic products by enabling purchase from a single supplier. Please ask your local sales representative to demonstrate how DYONICS ARTHROPAK Customized Procedure Kits can benefit your practice.

# Notes