

smith&nephew


FIRSTPASS◇


Suture Passer

Designed for consistent and secure suture passing



FIRSTPASS◇ Suture Passer






Key features

Pass and retrieve suture simultaneously

- Optional trap eliminates suture shuttling step
- Round, pre-curved nitinol needle provides added strength and reliable targeting
- Suture is carried in a recess on the needle which helps maintain suture integrity

Deep bite


- Jaw opens up to 10mm and has an 18mm bite depth




Round pre-curved needle for reliability and strength



Recess at the needle tip protects and carries the suture



Trap in top jaw captures suture loop



FirstPass can be used with or without the suture trap

Ordering information

FIRSTPASS Suture Passer

Reference #	Description
22-4035	FIRSTPASS SUTURE PASSER
22-4036	FIRSTPASS NEEDLE AND SUTURE CAPTURE, 5 PER BOX

ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735
USA

www.smith-nephew.com

Order Entry: 1-800-343-5717
Order Entry Fax: 1-888-994-2782

© 2015 Smith & Nephew, Inc.
◇Trademark of Smith & Nephew. Reg. US Pat. & TM Office.

P/N AT50-4008 Rev. C 02/15