smith&nephew

PERI-LOC◇ VLP

Variable-Angle
Locked Plating System

Distal Fibula Variable-Angle Locking Plates



PERI-LOC◇ VLP Locked Plating System

Indications

The PERI-LOC VLP Distal Fibula Plates are indicated for the treatment of fractures of the distal fibula.


The PERI-LOC VLP 3.5mm One-Third Tubular Locking Plates are indicated for the treatment of fractures, non-unions, and osteotomies of the medial malleolus, fibula, distal ulna, olecranon, calcaneus and metatarsals.


Contraindications

The PERI-LOC VLP Distal Fibula Plates are contraindicated for fractures with extreme metaphyseal comminution or dissociation of the articular segment from the shaft.


Design Features and Benefits

Polyaxial Locking Plate
Each PERI-LOC VLP screw hole contains five separate tabs that engage with the threads of the locking screw head to form a fixed angle construct. Locking screws can be angled and locked up to 15° in any direction allowing for the creation of customized, multi-directional locked plating constructs. Each screw hole accepts 3.5mm Cortex, 3.5mm Locking and/or 5.0mm Osteopenia Screws.






Low Profile Implants
The profile of the PERI-LOC VLP locking plate is designed to ensure low profile fracture fixation in areas of minimal soft tissue coverage such as periarticular zones. All screws in the PERI-LOC VLP Locked Plating System have a low head profile to further reduce the potential for soft tissue irritation in these sensitive areas.



Optimal Plate Contour
The optimized contour of the PERI-LOC® VLP locking plate facilitates both fracture reduction and stabilization as it is compressed to bone. This "compression contouring" feature is essential to the implant's ability to resist torque and bending during fracture healing. Once securely fixed in place, the plate produces a constant buttress effect to the fracture site to prevent loss of reduction and to enhance overall fracture fixation.



Implant Features

3.5mm Lateral Distal Fibula Locking Plate

- Distal screw cluster and low plate profile provide stable periarticular fixation
- 2.0mm – 1.7mm proximal to distal plate thickness transition
- Left and right specific



3.5mm PosteroLateral Distal Fibula Locking Plate

- Scalloped for syndesmotic screw placement outside the plate without compromising plate position
- 1.5mm – 0.9mm proximal to distal plate thickness transition
- Left and right specific
- 8° distal helical twist to accommodate the posterolateral anatomy of the distal fibula
- Rounded distal edge to minimize peroneal nerve irritation



3.5mm One-Third Locking Tubular Plate

- Low-profile buttress plate for fractures of the distal fibula
- Consistent 1.5mm plate thickness



Catalog Information

Cat. No.	Description
7282-0805	PERI-LOC® VLP 3.5mm Posterolateral Distal Fibula Locking Plate, 5 Hole Left 62mm
7282-0806	PERI-LOC VLP 3.5mm Posterolateral Distal Fibula Locking Plate, 6 Hole Left 74mm
7282-0807	PERI-LOC VLP 3.5mm Posterolateral Distal Fibula Locking Plate, 7 Hole Left 86mm
7282-0905	PERI-LOC VLP 3.5mm Posterolateral Distal Fibula Locking Plate, 5 Hole Right 62mm
7282-0906	PERI-LOC VLP 3.5mm Posterolateral Distal Fibula Locking Plate, 6 Hole Right 74mm
7282-0907	PERI-LOC VLP 3.5mm Posterolateral Distal Fibula Locking Plate, 7 Hole Right 86mm
7282-1003	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 3 Hole Left 59mm
7282-1004	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 4 Hole Left 71mm
7282-1005	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 5 Hole Left 83mm
7282-1007	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 7 Hole Left 107mm
7282-1009	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 9 Hole Left 131mm
7282-1011	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 11 Hole Left 155mm
7282-2003	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 3 Hole Right 59mm
7282-2004	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 4 Hole Right 71mm
7282-2005	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 5 Hole Right 83mm
7282-2007	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 7 Hole Right 107mm
7282-2009	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 9 Hole Right 131mm
7282-2011	PERI-LOC VLP 3.5mm Lateral Distal Fibula Locking Plate 11 Hole Right 155mm
7282-3005	PERI-LOC VLP 3.5mm One-Third Tubular Locking Plate 5 Hole 62mm
7282-3006	PERI-LOC VLP 3.5mm One-Third Tubular Locking Plate 6 Hole 74mm
7282-3007	PERI-LOC VLP 3.5mm One-Third Tubular Locking Plate 7 Hole 86mm
7280-3008*	PERI-LOC VLP 3.5mm One-Third Tubular Locking Plate 8 Hole 98mm
7280-3010*	PERI-LOC VLP 3.5mm One-Third Tubular Locking Plate 10 hole 122mm
7280-3012*	PERI-LOC VLP 3.5mm One-Third Tubular Locking Plate 12 Hole 146mm

* Available sterile only

Orthopaedic Trauma &
Clinical Therapies
Smith & Nephew, Inc.
1450 Brooks Road
Memphis, TN 38116
USA

Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders/Inquiries: 1-800-238-7538

www.smith-nephew.com

®Trademark of Smith & Nephew. Reg. US Pat. & TM Off.

30023405028a 7118-1194 11/07